Hazardous Waste Disposal Container Guide

Chemical Waste 1 gallon	Chemical Waste 5 gallons	Chemical Sharps Waste 2 gallons	Biohazardous Waste 20 gallons	Pathology Waste 44 gallons	Biohazardous Sharps Waste 2 gallons	I	
			BIOHAZARD WASTE	PATHOLOGY WASTE INCINERATION ONLY			
		Use For:					
Liquid or Solid Chemical Waste Use with secondary containment	Solid Chemical Waste	Chemically Contaminated Sharps	Solid Biohazardous Waste	Pathology Waste	Biohazardous Contaminated Sharps	E	
	•	-	Exam	ples:	*		
Methanol, formaldehyde, solid chemical powders, silica gel, etc.	Contaminated debris, gloves, bottles, glassware, Eppendorf tubes, ethidium bromide	Contaminated hypodermic needles, razors, scalpel blades, etc.	Infectious agents, bacteria, viruses, spores, etc.	Carcasses, recognizable body parts, transgenic animals, etc.	Contaminated or non-contaminated hypodermic needles, razors, scalpel blades, etc.	b	
	gels, etc.					1	
	gels, etc.		Dispose	within:			
6 months	gels, etc. 6 months	6 months	Dispose 7 days	within: 7 days	7 days of being 2/3 full		
6 months		6 months	•	7 days	7 days of being 2/3 full		



Clean **Broken Glass** 5 gallons

Locations

Scrap Metal Shavings and **Powder Waste** 5 gallons

